

# PROPEX<sup>®</sup> II

NEW GENERATION APEX LOCATOR WITH FULL COLORED DISPLAY



SHARPEN YOUR ROOT CANAL VISIBILITY

**DENTSPLY**  
MAILLEFER

# PROPEX<sup>®</sup> II




REF A 1028

## PRODUCT FEATURES

- Large colored display to facilitate the file tracking
- Progressive sound control for dual control
- Powered by rechargeable battery
- Latest multi-frequency technology enables location of the apex in most types of canal conditions
- Fully automatic device


## UPPER CANAL ZONE

The cursor on the tooth icon indicates the progression of the file in the canal




## APICAL ZONE


Extended apical zone from 0,9 to 0,0



When apex is reached a solid tone is emitted



When the cursor on the tooth icon shows over this indicates that the file has passed the apex



0459 Dentsply Maillefer • CH-1338 Ballaigues • Switzerland • [www.dentsplymailefer.com](http://www.dentsplymailefer.com)



ENDO4YOU

XX