

Do it yourself...

Ufi Gel hard Cartridge

first time in cartridge

Cartridge

- Art. No. 2215 Set 80 g cartridge, 10 ml adhesive, mixing tips type 8
- Art. No. 2218 Intro set 80 g cartridge, Dispenser type 2, 10 ml adhesive, mixing tips type 8
- Art. No. 2216 Refill pack 80 g cartridge, mixing tips type 8
- Art. No. 2217 Refill pack 10 ml adhesive



Handmix

- Art. No. 2210 Set 60 g powder, 40 ml liquid, 20 ml conditioner, accessories
- Art. No. 2211 Refill pack 60 g powder
- Art. No. 2212 Refill pack 40 ml liquid
- Art. No. 2213 Refill pack 20 ml conditioner



Available from:



VC 8400 2215 GB 027

direct chairside hard relining

now as simple as an impression

VOCO

P.O. Box 767 · 27457 Cuxhaven · Germany · Tel. +49 (4721) 719 - 0 · Fax +49 (4721) 719 - 140 · www.voco.com

creative in research

VOCO

creative in research



Ufi Gel hard Cartridge

Ufi Gel hard C is a PMMA-like material for permanent direct chairside hard relinings of full and partial dentures. It also allows lengthening of denture rims and adaptation of dentures layer by layer over the years. Relinings with this new material do not differ from laboratory processed pieces in their

hardness, finish or cleaning. The accurate direct impression with **Ufi Gel hard C** avoids loss of detail and the subsequent time-consuming corrections.

...more economical for you


- complete relining chairside in one session
- simple, time and material-saving application from the cartridge
- less organisational hassle in the clinic
- no costs for impression material or laboratory work
- better customer service (one short session, no waiting time without denture)
- easy and fast to polish

...more comfortable for your patient


- biocompatible (methylmethacrylate-free), neutral taste and odour
- low heat development (~37°C)
- bubble-free application, ensures a smooth surface
- durable tight bond to the denture material
- excellent flow properties for a precise fit
- aesthetic and stable colour



... as simple as an impression



Polymerisation temperature




Source: VOOCO in-house test, 1999
1 ml relining material after application to denture material

Colour change after 1 week

Before			
Storage in tea			
Storage in coffee			
	Ufi Gel hard C	A	B

hard chairside relining in only 10 minutes...



hard chairside relining: a win-win therapy